


# Comprehensive Chemokine Services

- One of the most comprehensive ranges of human and mouse chemokines in the native and site-specifically labelled form
- High purity, endotoxin free products
- Custom synthesis of human and non-human wild-type or labelled
   shame kines for products not available in the catalogue.





### Comprehensive Chemokine services

One of the most comprehensive ranges of human and mouse chemokines in the native and site- specifically labelled form

Proprietary skills in long chain peptide synthesis produce high purity, endotoxin free products. The selective manufacturing process is highly amenable to the manufacture of labelled products, ensuring that not only is a single label incorporated, but also at a specific site in the sequence. We provide site-specifically labelled chemokines and complement proteins with biotin and AlexaFluor647®1 from stock.

#### Our AlexaFluor®1 chemokines offer the following advantages:

- A powerful and unique reagent which absorbs at 650 nm and emits at 665 nm
- The highly selective nature of the synthesis process ensures only labelling at a single site in the protein sequence
- Wide ranging applications e.g. receptor binding studies<sup>2</sup> and flow cytometry<sup>3</sup>

### Chemokines

Systematic Name	Chemokine	HUMAN CHEMOKINES Chemokine Receptors	Native	Biotin	AlexaFluor647®1	
CCL1	h-I-309	CCR8	CN-7	CB-7	CAF-7	
CCL2	h-MCP-1	CCR2	CN-2	CB-2	CAF-2	
CCL3	h-MIP-1 alpha	CCR1/CCR5	CN-1	CB-1	CAF-1	
CCL4	h-MIP-1 beta	CCR5	CN-23	CB-23	CAF-23	
CCL5	h-RANTES	CCR1/CCR3/CCR4/CCR5	CN-8	CB-8	CAF-8	
CCL8	h-MCP-2	CCR1/CCR2/CCR3/CCR5	CN-18	CB-18	CAF-18	
CCL11	h-Eotaxin	CCR3	CN-3	CB-3	CAF-3	
CCL15	h-MIP-5	CCR1/CCR3	CN-19	CB-19	CAF-19	
CCL17	h-TARC	CCR4	CN-16	CB-16	CAF-16	
CCL18	h-PARC	unknown	CN-21	CB-21	CAF-21	
CCL19	h-MIP-3 beta	CCR7	CN-6	CB-6	CAF-6	
CCL20	h-MIP-3 alpha	CCR6	CN-5	CB-5	CAF-5	
CCL22	h-MDC	CCR4	CN-4	CB-4	CAF-4	
CCL25	h-TECK	CCR9	CN-15	CB-15	CAF-15	
CCL28	h-MEC	CCR10	CN-20	CB-20	CAF-20	
CXCL8	h-IL-8	CXCR1/CXCR2	CN-9	CB-9	CAF-9	
CXCL10	h-IP-10	CXCR3	CN-10	CB-10	CAF-10	
CXCL11	h-ITAC	CXCR3	CN-13	CB-13	CAF-13	
CXCL12	h-SDF-1 alpha	CXCR4	CN-11	CB-11	CAF-11	
CXCL14	h-BRAK	unknown	CN-17	CB-17	CAF-17	
CX3CL1	h-Fractalkiine	CX3CR1	CN-14	CB-14	CAF-14	
MOUSE CHEMOKINES						
Systematic Name	Chemokine	Chemokine Receptors	Native	Biotin	AlexaFluor647®1	
mCCL2	mouse MCP-1	CCR2	CN-52	CB-52	CAF-52	
mCXCI 12	mouse SDF-1 alpha	CXCR4	CN-50	CB-50	CAF-50	

CX3CR1

COMPLEMENT PROTEINS							
Complement Protein	Native	Biotin	AlexaFluor647®1				
hC3a	CN-91	CB-91	CAF-91				
hC5a	CN-90	CB-90	CAF-90				

CN-51

CB-51

CAF-51

Native available in 20µg, 100µg, 1mg Biotin and Alexa647®1 available in 10µg, 100µg

mCX3CL1

Email: chemokines@almacgroup.com

mouse Fractalkine

<sup>&</sup>lt;sup>1</sup> AlexaFluor 647®- labelled chemokines are sold under license for research use only patents or patent applications owned by Molecular Probes


<sup>&</sup>lt;sup>2</sup> Strong et al (2006). "Synthetic chemokines directly labelled with a fluorescent dye as tools for studying chemokine and chemokine receptor interactions". Eur.Cytokine Netw., 17 (1), 49-59

<sup>&</sup>lt;sup>3</sup> Le Brocq et al (2014). "Chemokines as novel and versatile reagents for flow cytometry and cell sorting".

J Immunol, 192, 6120-6130

## Complementary peptide synthesis services:

- · Custom Synthesis
- Challenging Targets
- Protein Engineering
- · Peptide Conjugation
- · Peptide Labelling
- Analytical Support.



### Get in touch

### UK

Almac Group 5 The Fleming Building Edinburgh Technopole Milton Bridge Nr Penicuik EH26 0BE United Kingdom

chemokines@almacgroup.com +44 28 3839 5717

#### US

Almac Group (US Headquarters) 25 Fretz Road Souderton, PA 18964 United States of America

chemokines@almacgroup.com +1 215 660 8500