

**h-IP-10 (AF647®), synthetic****100 µg**

Catalogue number: CAF-10

**Almac Peptide and Protein Technologies****Chemokines**

Custom Peptides

Site-Specific protein labelling

Modified Histones

Ubiquitylated peptides

**Product Information**

<b>Product Name:</b>	h-IP-10 (AF647®) CXCL10
<b>Synonyms:</b>	CXCL10 (Alexa 647®) human, 10 kDa interferon-gamma-induced protein, Gamma-IP10, crg-2, Gamma-Interferon Inducible Protein 10.
<b>Swiss Prot Accession No:</b>	P02778
<b>Sequence:</b>	VPLSRTVRCT CISISNQPVN PRSLEKLEII PASQFCPRVE IIATMKKKGE KRCLNPESKA IKNLLKAVSK ERSKRSP-Linker- K(AF647®)
<b>Modifications:</b>	AF647® incorporated at indicated position in a molar ratio of dye:peptide of 1:1. Linker used is 8-amino-3,6-dioxaoctanoic acid.

**Specification**

<b>HPLC Purity:</b>	>95% area/area at 214nm
<b>Determined Mass:</b>	9.7 kDa
<b>Formulation / Appearance:</b>	Blue powder lyophilized from water/acetonitrile/ TFA to generate the TFA salt form.  Peptide content determined by UV analysis at 280nm.

**Preparation and Storage**

<b>Reconstitution / Storage:</b>	It is recommended that unopened vials are stored at -20 °C to -70 °C for periods of up to 18 months. Avoid repeat freeze-thaw cycles.  It is recommended vials be centrifuged prior to opening.  Water can be used to prepare stock solutions of 20 µmol.L-1. Stock solutions with up to 30% DMSO/water can also be prepared.
----------------------------------	---

*Not fully tested. For research use only. Not for use in human diagnostic or therapeutic procedures.*