

T+44 (0)28 3839 5717 • E almacshop@almacgroup.com • W shop.almacgroup.com

h-TECK (Biotin), synthetic

10 µg

Catalogue number: CB-15

Almac Peptide and Protein Technologies

Chemokines

Custom Peptides

Site-Specific protein labelling

Modified Histones

Ubiquitylated peptides

Product Information

Product h-TECK (Biotin)

Name: CCL25

Synonyms: CCL25 (Biotin) human, Small inducible cytokine A25

(Precursor) truncate.

Swiss Prot

Accession No:

Sequence: Pyr-GVFEDCCLA YHYPIGWAVL RRAWTYRIQE

VSGSCNLPAA IFYLPKRHRK VCGNPKSREV QRAMKLLDAR

NK(Linker-Biotin)VF

Modifications: Biotin incorporated at indicated position in a molar ratio of

biotin:peptide of 1:1. Linker used is 8-amino-3,6-dioxaoctanoic

acid.

O15444

Specification

HPLC Purity: >95% area/area at 214nm

Determined

Mass:

8.9 kDa


Formulation /

White powder lyophilized from water/acetonitrile/ TFA to

Appearance: generate the TFA salt form.

Peptide content determined by UV analysis at 280nm.

Preparation and Storage



Reconstitution / Storage:

It is recommended that unopened vials are stored at -20 °C to -70 °C for periods of up to 18 months. Avoid repeat freeze-thaw

cycles.

It is recommended vials be centrifuged prior to opening.

Water can be used to prepare stock solutions of 20 μ mol.L⁻¹. Stock solutions with up to 30% DMSO/water can also be

prepared.

Not fully tested. For research use only. Not for use in human diagnostic or therapeutic procedures.